

AMXD Advanced Mission Extender Device

Component Data Sheet

Omni Part Number: **CB-10 & CB-12 & CB-16**

Collection Bag 800ML Capacity (CB-10)
Collection Bag 1200ML Capacity (CB-12)
Collection Bag 1600ML Capacity (CB-16)



Disposable

Gels on contact with liquid

One way valve

Size of 800ML Capacity (CB-10): 11.5" L x 5" W.

Contents: Super Absorbent Polymer sheets

Environmentally Friendly

Safe, Non-toxic



Packaged

Omni Measurement System, Inc.

115 Catamount Drive, Milton, VT 05468 802-891-5500 Sales@OmniAMXD.COM