

AMXD Advanced Mission Extender Device

Component Data Sheet

Omni Part Number: **MG-1-B & MG-2-B**

Male Undergarments Size 1 & Size 2 (BOXER STYLE)



Environmentally Sealed



Reusable

Sizing:

Choose Size 1 (for waist size between 30" – 34")

Choose Size 2 (for waist size between 36" – 40")

Washable

Style: Boxer

100% Cotton

Omni Measurement System, Inc.

115 Catamount Drive, Milton, VT 05468 802-891-5500 Sales@OmniAMXD.COM