

AMXDmax

Removeable Hose



RHA-21B-2-X

Standard 12", 18", or 24"



ESC Compatible



Features and Benefits

- Flame Retardent Nomex covering
- Kink resistant corrugated tubing
- Separates from Male Cup / Female Pad for easy cleaning and storage
- Sensors in the hose detect urine for fully automated bladder relief
- Available in 18", 24", or ESC Compatible Versions
- Certified Safe to Fly

Catalog Number

RHA-21B-2-18* - 18" Length

RHA-21B-2-24* - 24" Length

RHA-21B-2-12* - ESC Compatible

GSA

SCU-ABP-B (18" Hose)

SCU-ABP-B (24" Hose)

SCU-ABP-B (12" Hose)

Description

The removeable hose is the detachable hose for use with the Male Cup or Female Pad. It contains the sensors to detect urine when attached to the Male Cup / Female PAD.

Related Products

Male Cup w/ Removeable Hose, Female Pad w/ Removeable Hose, Through-Suit Connector, Collection Bag

Specifications

Length*	18" / 24" / 12" (ESC)
Weight	1.22oz / 1.5oz
Operating Temp.	14° F to 122° F -10° C to 50° C
Storage Volume	6" x 9" x 3" 152mm x 229mm x 76mm
Storage Temp.	-4° F to 158° F -20° C to 70° C

* specify Male/Female when ordering



Omni Medical Systems, Inc.

808 Hercules Drive, Colchester, VT 05446

☎: (802)497-2253 ☎: (802)497-3601

www.amxdmax.com