

# **AMXD    Advanced Mission Extender Device**

## **Component Data Sheet**

Omni Part Number:    **SCU-10 & SCU-10A**

**System Control Unit**



The system control unit is utilized with the IFP/IMC and Collection Bag to provide aircrew with in-flight bladder relief.

The SCU-10 includes:

-One Pump unit containing air inflation systems and liquid pumping hardware.

Dimensions :

65mm x 44.5mm x 145mm

Weight :

.756 lbs without battery pack loaded

**Omni Measurement System, Inc.**

115 Catamount Drive, Milton, VT 05468 802-891-5500 Sales@OmniAMXD.COM