

AMXDmax

Rechargeable Battery Pack



SCU-RBP-2



Features and Benefits

- Extended capacity lithium polymer rechargeable battery pack
- Capacity standard 420 and 630 for 50% greater capacity
- Uses standard AMXDmax battery charger
- Certified Safe to Fly

Catalog Number

SCU-RBP-2 420*
SCU-RBP-2 630*

GSA

SCU-RBP (-2 Specify 420)
SCU-RBP (-2 Specify 630)

Description

The Rechargeable Battery pack is the battery used with the AMXDmax Control Unit. It is available in 2 sizes.

Related Products

Disposable Pack (Holder)

Specifications

Length	2.38" (60mm)
Width	1.2" (30mm)
Height	0.55" (14mm) / 0.71" (18mm)
Weight	0.9oz / 1.3oz
Operating Temp.	14° F to 122° F -10° C to 50° C
Storage Volume	2.38" x 1.2" x 0.55" / 0.71" 60mm x 30mm x 14mm / 18mm
Storage Temp.	-4° F to 158° F -20° C to 70° C

* specify when ordering



Omni Medical Systems, Inc.

808 Hercules Drive, Colchester, VT 05446

☎: (802)497-2253 ☎: (802)497-3601

www.amxdmax.com